

Martin Lindquist

Department of Biostatistics
Johns Hopkins
Bloomberg School of Public Health

Tor Wager

Department of Psychology and
Neuroscience and the
Institute for Cognitive Science
University of Colorado, Boulder

Image Formation

Signal Formation

- The subject is placed into the MR scanner.
 - Nuclei of ^1H atoms align with the magnetic field.
 - The nuclei precess about the field at similar frequencies, but at a random phase.
 - Net longitudinal magnetization in the direction of field.
- Within a slice, a radio frequency (RF) pulse is used to align the phase and ‘tip over’ the nuclei.
 - Causes the longitudinal magnetization to decrease, and establishes a new transversal magnetization.

Signal Formation

- After the RF pulse is removed, the system seeks to return to equilibrium.
 - The transverse magnetization disappears ([transversal relaxation](#)), and the longitudinal magnetization grows back to its original size ([longitudinal relaxation](#)).
 - Longitudinal relaxation: exponential growth described by time constant T1.
 - Transverse relaxation: exponential decay described by time constant T2.
- During this process a signal is created that can be measured using a receiver coil.

Slice Selection

- Most structural MRI and fMRI scans involve the construction of a three dimensional image from a set of two-dimensional slices.






Image Formation

- Imagine a brain slice split into a number of equally sized **volume elements** or **voxels**.



$\rho(x,y)$



Gradients

- The measured signal combines information from the whole brain:

$$S(t) = \iint \rho(x, y) dx dy$$

- A **magnetic field gradient** is used to sequentially control the spatial inhomogeneity of the magnetic field, so each measurement can be expressed:

$$S(k_x, k_y) = \iint \rho(x, y) e^{-i2\pi(k_x x + k_y y)} dx dy$$

K-space


- The measurements are acquired in the frequency-domain (**k-space**).
- By making measurements for multiple values of (k_x , k_y) we can gain enough information to solve the inverse problem and reconstruct $\rho(x, y)$.
- We can use the **inverse Fourier transform (IFT)**:

$$\rho(x, y) = \iint S(k_x, k_y) e^{i2\pi(k_x x + k_y y)} dk_x dk_y$$




K-space Measurements

- In practice, data measurements are made *discretely* over a *finite region*.
 - Use discrete Fourier transforms.
- The number of k-space measurements we make influences the spatial resolution of the image.
 - Need enough measurements to solve inverse problem.





16 unknowns



4 unknowns

EPI and Spirals



Magnitude Images




- The measured k-space data is complex valued.
 - Hence, the measurement at each voxel is complex.
- We typically work with magnitude images, or

$$|\rho(x, y)| = \sqrt{\rho_R(x, y)^2 + \rho_I(x, y)^2}$$

where ρ_R and ρ_I are the real and imaginary parts of the k-space measurement.



Image Formation



End of Module



@fMRIstats