

Martin Lindquist

Department of Biostatistics
Johns Hopkins
Bloomberg School of Public Health

Tor Wager

Department of Psychology and
Neuroscience and the
Institute for Cognitive Science
University of Colorado, Boulder

fMRI Artifacts & Noise

fMRI Noise

- BOLD fMRI signal contains multiple sources of noise related to the hardware and the participants themselves.
- Sources of noise:
 - Thermal motion of free electrons in the system
 - Gradient and magnetic field instability
 - Head movement and its interactions with magnetic field
 - Physiological effects, including heartbeat and respiration
- Appears in data as:
 - High-frequency ‘spikes’
 - Image artifacts and distortions
 - Low-frequency (slow) drift and periodic fluctuation over time

Noise and artifact mitigation


- Acquisition:
 - Good quality control processes for scanner
 - Choose acquisition sequence/parameters for your goal
 - Specialized sequences (e.g., spin-echo, simultaneous multi-slice, z-shimming).
 - Minimize head movement

- Analysis:
 - Look at data as it is acquired to check for issues
 - Outlier / artifact identification and correction
 - Preprocessing to adjust for head movement, drift
 - Helpful statistical procedures (e.g., robust regression, hierarchical modeling)
 - Low-frequency (slow) drift and periodic fluctuation over time


Acquisition artifacts: Look at your data!

- Five issues you should check for and avoid


1. Coverage not what you intended
2. RF noise and malformed images
3. Transient gradient artifacts/spikes
4. Ghosting
5. Susceptibility artifacts/dropout
6. Task-correlated head movement




Task-correlated movement




5




1



2




3



4

Ghosting image courtesy of Luis Hernandez-Garcia


Acquisition artifacts: Look at your data!



- Five issues you should check for and avoid
 1. Coverage not what you intended
 2. RF noise and malformed images
 3. Transient gradient artifacts/spikes
 4. Ghosting
 5. Susceptibility artifacts/dropout
 6. Task-correlated head movement
 - All fMRI data contain some artifacts
 - It is very difficult to deal with bad artifacts during analysis
 - Try to avoid them during acquisition


Drift

- Slow changes in voxel intensity over time (low-frequency noise) is present in the fMRI signal.
- Scanner instabilities are a main cause of drift, as drift has been seen in cadavers, but aliased physiological noise is also important.
- We need to account for drift both in preprocessing and when we conduct statistical analyses.



Issues

- Drift can have serious consequences:
 - Experimental conditions that vary slowly may be confused with drift.
 - Experimental designs should use high frequencies (more rapid alternations of stimulus on/off states).
- Bad Design:



Motion

- Subject motion during the experiment can also give rise to serious problems.
- Typically motion correction is performed in the pre-processing stages of the analysis.
- However, ‘spin-history’ artifacts may remain that cannot be removed.
 - This is caused in part by through-plane motion and complex interactions with the magnetic field
 - Attempts to account for it are often made in during statistical analysis, but these are imperfect.




Physiological Noise

- Respiration and heart beat give rise to noise at particular frequencies.
- It can potentially be modeled in statistical analysis, but if the TR is too low there will be problems with *aliasing*.
- For standard TR values (~ 2s) this type of noise is difficult to remove and is often left in the data giving rise to temporal autocorrelations.




Aliasing


- Periodic signals that occur more rapidly than the sampling rate will often be *aliased* back to lower frequencies.



Signal (blue) sampled at near its periodicity. Observed signal is aliased, close to the ‘fundamental frequency’ (flat line).



Another example: Sampling at about $\frac{1}{2}$ the original frequency results in aliased periodicity about $\frac{1}{4}$ the original frequency.




Signals faster than $\frac{1}{2}$ the sampling rate—called the Nyquist frequency—will be aliased. To avoid aliasing, you must sample at least twice as fast as the fastest frequency in your signal.

Modeling fMRI Noise

- Some noise components can be removed prior to or during analysis, including low-frequency drift and images identified as outliers.
- However, it is impossible to remove or model all sources of noise, and therefore significant autocorrelation (non-independence of the observations over time) is usually present in the signal.
- In fMRI we typically use autoregressive (AR) or autoregressive moving-average (ARMA) processes to model the autocorrelation.

Spatio-temporal Behavior

- The spatiotemporal behavior of these noise processes is complex.



Spatial maps of the model parameters from an AR(2) model estimated for each voxel's noise data.

End of Module



@fMRIstats