

Martin Lindquist

Department of Biostatistics
Johns Hopkins
Bloomberg School of Public Health

Tor Wager

Department of Psychology and
Neuroscience and the
Institute for Cognitive Science
University of Colorado, Boulder

Basic MR Physics

Magnetic Resonance Imaging



An MR scanner consists of an electromagnet with a very strong magnetic field (1.5 - 7.0 Tesla)

Earth's magnetic field = 0.00005 Tesla

3 Tesla is 60,000 times stronger than the Earth's magnetic field.


What MRI Measures

- MRI is an extremely versatile imaging modality that can be used to study both brain structure and brain function.
- Both structural and functional MRI images are acquired using the same scanner.
- Different types of brain images can be generated to emphasize **contrast** related to different tissue characteristics.


Principles of MRI

- All magnetic resonance imaging techniques rely on a core set of physical principles.
- To understand we must begin by studying a **single atomic nuclei** and illustrate its impact on the generated MR signal.
- In particular we focus on **hydrogen** atoms consisting of a single proton.

Protons can be viewed as positively charged spheres which are always spinning. They give rise to a net magnetic moment along the axis of the spins.




Protons can be viewed as positively charged spheres which are always spinning. They give rise to a net magnetic moment along the axis of the spins.




Net Magnetization


- We cannot measure the magnetization of a single proton using MR, instead we measure the **net magnetization** of all nuclei within a volume.
- The net magnetization M can be viewed as a vector with two components.
 - A **longitudinal component** parallel to the magnetic field.
 - A **transverse component** perpendicular to the field.









In the absence of an external magnetic field, the nuclear magnetic moments are randomly oriented. There is no net magnetization.




When placed in a strong magnetic field, the nuclei align with the field. This creates a net longitudinal magnetization in the direction of the field.




The nuclei precess about the field with an angular frequency determined by the [Larmor frequency](#) but at a random phase.



A radio frequency (RF) pulse is used to align the phase and ‘tip over’ the nuclei. This causes the longitudinal magnetization to decrease, and establishes a new transversal magnetization.



A radio frequency (RF) pulse is used to align the phase and ‘tip over’ the nuclei. This causes the longitudinal magnetization to decrease, and establishes a new transversal magnetization.



Relaxation


- After the RF pulse is removed, the system seeks to return to equilibrium.
- The transverse magnetization starts to disappear (transversal relaxation), and the longitudinal magnetization grows back to its original size (longitudinal relaxation).
- During this process a signal is created that can be measured using a receiver coil.

Relaxation

- **Longitudinal Relaxation** is the restoration of net magnetization along the longitudinal direction as spins return to their parallel state.
 - Exponential growth described by time constant T1
- **Transverse Relaxation** is the loss of net magnetization in the transverse plane due to loss of phase coherence.
 - Exponential decay described by time constant T2




Longitudinal Relaxation Time




The restoration of longitudinal magnetization is described by a time constant T_1 .

Longitudinal Relaxation Time



Transverse Relaxation Time



The decay of magnetization due to interaction between nuclei is described by a time constant T_2 .

Image Contrast

- By altering how often we excite the nuclei (TR) and how soon after excitation we begin data collection (TE) we can control which characteristic is emphasized.
- The measured signal is approximately

$$M_0(1 - e^{-TR/T_1})e^{-TE/T_2}$$

where T_1 and T_2 are tissue properties.




Image Formation

- The goal of MRI is to construct an **image**, or a matrix of numbers that correspond to spatial locations.
- The image depicts the spatial distribution of some property of the nuclei within the sample.
- This could be the **density** of nuclei or the **relaxation time** of the tissues in which they reside.



Image Contrast



$$M_0(1 - e^{-TR/T_1})e^{-TE/T_2}$$

TE (echo time) -
the time between
excitation and data
collection.

fMRI Contrast

- T_2^* is the combined effect of T_2 and local inhomogeneities in the magnetic field.
- The scanner can be programmed to eliminate the effects of these inhomogeneities, or alternatively emphasize them.
- The latter types of procedures form the basis of BOLD fMRI.



Image Contrast

- Images can be produced that are sensitive primarily to T_1 , T_2 , or T_2^* .
- Because T_1 and T_2 vary with tissue type, they are able to represent boundaries between CSF, gray and white matter.
- Because T_2^* is sensitive to flow and oxygenation, it can be used to image brain function.

End of Module



@fMRIstats